

RK-05



Zna?ka: Dino-Lite
Kód produktu: RK-05

Short Description

Portable rack
Tilt and stay feature (30 degrees)
Fast assembly

Popis

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05 stand provides a multi-angle view with its 30 degrees inclinable pole. RK-05 is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05F with fine focus adjustment](#)

Specifikace

| | |
|------------------|---|
| Compatibility | |
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing | |

| | |
|----------------------|----------------------------------|
| Dimensions | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Weight | 611 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | Stand, HD-P1 holder |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €150,00 - €250,00 |

Product Gallery

