

MS15X

Brand: Dino-Lite
Product Code: MS15X



Short Description

High precision X/Y table

Description

The MS15X is an advanced XY table that can be used independently or in combination with the RK-06A or RK-10A column stands.

The MS15X X/Y table base with detachable rotating table is designed for use with all sorts of applications that require fine examination in an X- and Y-axis movement. With the 360 degree rotating table, one can easily carry out all kinds of precision examinations. The MS15X can be used stand-alone or attached to the RK-06A and RK-10A stands (with the included C-Clamps)

Specification

| | |
|------------------|---|
| Compatibility | |
| Compatible with: | All models; can be mounted on RK-06A/RK-10A |
| Housing | |
| Dimensions | 12.7cm (W) x 15.2cm (L) x 7.6cm (H) |
| Weight | 907 g |
| Features | |
| ESD safe | No |
| Information | |

| | |
|----------------------|---------------------------|
| Package contents | Base, stage, clamps |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €150,00 - €250,00 |

Product Gallery

