

RK-02

Marke: Dino-Lite

Artikelnummer: RK-02



Short Description

Gooseneck stand

Beschreibung

The RK-02 is a durable gooseneck stand solution designed for flexible and versatile usage with a Dino-Lite digital microscope. The C-clamp of the RK-02 fits to virtually any surface of up to 7 cm (~3 inch) thick, and can be mounted horizontally or vertically. With this sturdy stand, you can easily position the Dino-Lite by bending the 60 cm long arm and adjusting the rotating head. The base of the RK-02 can also be rotated providing even more freedom and it can be fixed with the quick release fastener for use on uneven surfaces. The RK-02 gooseneck arm is the ideal choice to flexibly use any of the Dino-Lite models, except the high magnification models.

Spezifikation

| | |
|------------------|--|
| Compatibility | |
| Compatible with: | All Dino-Lite models, except High Magnification models |
| Housing | |
| Dimensions | Gooseneck 60cm (L) x 1.5cm (D) |
| Weight | 1.5 kg |
| Features | |
| ESD safe | No |

| Information | |
|----------------------|---------------------------|
| Package contents | Stand, HD-P1 holder |
| Warranty information | 2 years European warranty |
| Accessory range: | Professional stands |
| Price range | €100,00 - €200,00 |

Product Gallery

