

# CA1070

**Marca:** Dino-Lite

**Codice del prodotto:** CA1070



## Short Description

Sturdy watertight carrying case

## Descrizione

Sturdy, watertight case for Dino-Lite microscopes. This heavy duty carrying case with foam inlay protects your Dino-Lite against dust, scratches and moisture. Watertight according to IP67 standard.

## Specifiche

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	Exterior: 210 x 167 x 90mm Interior: 186 x 123 x 75mm
Weight	570 g
Features	
ESD safe	No
Information	
Package contents	Case, foam inlay
Warranty information	2 years European warranty
Accessory range:	Various
Price range	€0 - €50,00

**Product Gallery**

