

# EarScope Basic (MEDL3E)



**Merk:** Dino-Lite Medical

**Productcode:** MEDL3E

## Korte omschrijving

640x480 Resolution

Basic model for otoscopy

20-30x Magnification

Medical Device Class 1

## Omschrijving

The EarScope Basic (MEDL3E) is a compact and safe solution to look at the ear. The cap is detachable for hygienical reasons. The MEDL3E has a 640x480 resolution and magnification of 20 ~ 30x. This model is delivered with disposable specula of 3, 4 and 5 mm (10 each).



USB 2.0



Adjustable  
~20x-30x



VGA resolution  
640 x 480



4 White LEDs  
non switchable

## Specificatie

Verlichting	
Licht / LED-type	White
Aantal leds	4

LED aan / uit schakelbaar:	Yes
Infrarood filter	IR cut-filter >650 nm
Diffusor beschikbaar	No
Polarisator	No
<b>Optiek</b>	
Vergroting	20-30x
Type lens	Glass with anti-reflection coating
<b>Sensor</b>	
Sensortype	CMOS
Resolutie	VGA (640x480)
Maximale framesnelheid	30 fps
<b>Compatibiliteit</b>	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
<b>Behuizing</b>	
Materiaal van de behuizing	Composite/plastic housing
Vergrotingsslot	No
Afmetingen	10.5cm (L) x 3.2cm (D)
Gewicht	95 g
Kabellengte	1.8m
<b>Kenmerken</b>	
Meting	No
Kalibratie	No
Microtouch-sensor	No
<b>Informatie</b>	
Inhoud van de verpakking	Microscope, carry pouch, 10x 3 mm specula, 10x 4 mm specula, 10x 5 mm specula, software CD, user manual
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745

Prijsklasse	€200,00 - €300,00
-------------	-------------------

## **Product galerij**