

# RK-02

**Merk:** Dino-Lite

**Productcode:** RK-02



## Korte omschrijving

Gooseneck stand

## Omschrijving

The RK-02 is a durable gooseneck stand solution designed for flexible and versatile usage with a Dino-Lite digital microscope. The C-clamp of the RK-02 fits to virtually any surface of up to 7 cm (~3 inch) thick, and can be mounted horizontally or vertically. With this sturdy stand, you can easily position the Dino-Lite by bending the 60 cm long arm and adjusting the rotating head. The base of the RK-02 can also be rotated providing even more freedom and it can be fixed with the quick release fastener for use on uneven surfaces. The RK-02 gooseneck arm is the ideal choice to flexibly use any of the Dino-Lite models, except the high magnification models.

## Specificatie

Compatibiliteit	
Compatibel met:	All Dino-Lite models, except High Magnification models
Behuizing	
Afmetingen	Gooseneck 60cm (L) x 1.5cm (D)
Gewicht	1.5 kg
Kenmerken	
ESD veilig	No

Informatie	
Inhoud van de verpakking	Stand, HD-P1 holder
Garantie-informatie	2 years European warranty
accessoires:	Professional stands
Prijsklasse	€100,00 - €200,00

### Product galerij

