

F25MM

Merk: Dino-Lite
Productcode: F25MM



Korte omschrijving

C-mount adapter for rigid endoscope

Omschrijving

C-mount adapter for rigid endoscopes with a focal length of 25mm. This C-mount adapter connects to rigid, straight through endoscopes with max. 32mm ocular diameter.

Specificatie

Compatibiliteit	
Compatibel met:	AM4023CT, AM4025X, AM7025X
Behuizing	
Afmetingen	48mm (H) x 52mm (W)
Gewicht	120 g
Kenmerken	
ESD veilig	No
Informatie	
Inhoud van de verpakking	Adapter
Garantie-informatie	2 years European warranty
accessoires:	Various
Prijsklasse	€250,00 - €350,00

Product galerij

