

# WF-20

**Merk:**

**Productcode:** WF-20



## **Korte omschrijving**

Stable transmission  
Network connectivity  
Long battery life  
Adaptable design  
Agility

## **Omschrijving**

Coupling with Dino-Lite AF series, the WF-20 is a Wi-Fi streamer adding further versatility, flexibility, and agility to Dino-Lite microscope. Thanks to stable image transmission and long battery life (more than 2,5 hours), the WF-20 is well suited but not limited to field applications. The WF-20 can also be configured as a wireless router for staying connected with internet or a network during usage.

The WF-20 allows AF microscope models to be used wirelessly in connection with an iOS (iPhone/iPad) or Android app on any tablet, smartphone or even in the DinoCapture 2.0 software for Windows computers. The WiFi streamer is ideal for field work or presentations and can stream from one Dino-Lite to multiple devices. When using the iOS app or DinoCapture, users can perform measurements and take advantage of the automatic magnification reading (AMR) with compatible Dino-Lite AMR and measurement models. The WF-20 can be used anywhere to create a WiFi signal for streaming the Dino-Lite live images. The free DinoConnect App for WF-10 can be downloaded from the App Store or Google Play.

For (non "AF") USB microscope models see the [WF-10](#) for wireless options.

---

[Download the WF-20 instruction manual](#)



## Specificatie

Compatibiliteit	
Besturingssysteem	Windows 7, 8, 10 & 11, iOS 8.0 or later, Android 4.4 or later
Software	DinoCapture 2.0 (Windows), DinoConnect App
Compatibel met:	AF4x15 & AF7x15 series
Behuizing	
Afmetingen	58mm (L) x 36mm (W) x 32mm (H)
Gewicht	57 g
Kenmerken	
ESD veilig	No
Frequency:	2.412GHz~2.485GHz
Informatie	
Inhoud van de verpakking	WF-20 WiFi streamer, user manual, mini USB charging cable, integrated Li-Ion battery
Garantie-informatie	2 years European warranty
accessoires:	Wireless
Prijsklasse	€250,00 - €400,00

## Product galerij

