

# RK-10A

**Marka:** Dino-Lite  
**Ürün Kodu:** RK-10A

## Short Description

High precision stand with quick release

## Açıklama

The Dino-Lite RK-10A is a sturdy and stable high-end stand solution that can be used with all the professional Dino-Lite models. It is constructed of resilient stainless steel and lightweight aluminum and offers a very precise fine-focus adjustment as well as a quick release function.

The quick release button allows for fast and convenient vertical movement and the fine-focus adjustment knob is especially useful at high magnification. This model also includes an adjustable safety stop to prevent the microscope from touching the specimen. Its adjustable horizontal arm has a 15cm range with 360-degree rotation and the additional arm (RK-10-EX included) increases its effective range to observe larger objects. With these features, the RK-10A is the ideal accessory for precise orientation and positioning for maximum productivity and ease of use in a wide range of environments including industrial, scientific and laboratory settings.

## [AVAILABLE EXTENSIONS FOR THE RK-10A STAND](#)

---

### **RK-10-EX**



additional horizontal arm (included with RK-10A)

**RK-10-FX**




flexible arm extension

**RK-10-VX**



vertical arm extension

**RK-10-PX**



*(Temporarily unavailable)*

### **Özellikler**

|                  |                                |
|------------------|--------------------------------|
| Compatibility    |                                |
| Compatible with: | All Dino-Lite models           |
| Housing          |                                |
| Dimensions       | 15cm (W) x 22cm (L) x 23cm (H) |
| Weight           | 1.10 kg                        |
| Features         |                                |
| ESD safe         | Yes                            |

## Information

|                      |  |
|----------------------|--|
| Package contents     | Stand, HD-P1 holder, RK-10-EX horizontal arm |
| Warranty information | 2 years European warranty                    |
| Accessory range:     | Professional stands                          |
| Price range          | €300,00 - €450,00                            |

## Product Gallery



